# Tetric® N-Ceram Tetric® N-Flow Tetric® N-Bond N-Etch

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English

### Descriptions

- **N-Etch** is a phosphoric acid gel for enamel etching and dentin conditioning.
- **Tetric N-Bond** is a light-curing, single-component bonding agent for enamel and dentin bonding in combination with the total-etch technique.
- **Tetric N-Flow** is a flowable, light-curing, radiopaque nano-hybrid composite for restorative treatment and the cementation of translucent ceramic and composite restorations.
- **Tetric N-Ceram** is a light-curing, radiopaque nano-hybrid composite for restorative treatment.

Tetric N-Ceram, Tetric N-Flow and Tetric N-Bond cure with light in the wavelength range of 400–500 nm (blue light).

### Shades

**Tetric N-Ceram** is available in the following shades:

*Enamel shades:* A1, A2, A3, A3.5, A4, B2, B3, C2, C3, D3

*Dentin shades:* A3.5 Dentin, B2 Dentin

*Highly translucent Incisal shade:* T

*Bleach shades:* Bleach Light (L), Bleach Incisal (I), Bleach Medium (M)

**Tetric N-Flow** is available in the following shades:

*Enamel shades:* A1, A2, A3, A3.5, A4, B2

*Dentin shades:* A3.5 Dentin

*Highly translucent Incisal shade:* T

*Bleach shades:* Bleach Light (L), Bleach Incisal (I)

### Composition

- **N-Etch** contains phosphoric acid (37 wt. % in water), thickeners and pigments.
- **Tetric N-Bond** contains phosphoric acid acrylate, HEMA, Bis-GMA, urethane dimethacrylate, ethanol, film-forming agent, catalysts and stabilizers.
- **Tetric N-Flow** contains 36 wt.% dimethacrylates (including TEGDMA), 63 wt.% fillers (barium glass, ytterbium trifluoride, highly dispersed silica and mixed oxide) and 1 wt.% catalysts, stabilizers and pigments. The total content of inorganic fillers is 39 vol.%. The particle size of inorganic fillers is between 40 and 3000 nm.
- **Tetric N-Ceram** consists of dimethacrylates (19-20 wt.%). The fillers contain barium glass, ytterbium trifluoride, mixed oxide and copolymers (80-81 wt.%). Additives, catalysts, stabilizers and pigments are additional contents (< 1 wt.%). The total content of inorganic fillers is

55–57 vol.%. The particle size of inorganic fillers is between 40 nm and 3000 nm.

### Indications

- **N-Etch** is used in the enamel etch or total-etch-technique in conjunction with adhesives, composites, luting materials and fissure sealants.
- **Tetric N-Bond**
  - Adhesive for direct, light-curing and dual-curing composite and compomer restorations
  - Adhesive for indirect all-ceramic and composite restorations (inlays, onlays, veneers)
- **Tetric N-Flow**
  - As initial layer under Class I and II restorations
  - Anterior restorations (Classes III, IV)
  - Class V restorations (cervical caries, root erosion, wedge-shaped defects)
  - Small restorations of all types
  - Extended fissure sealing
  - Splinting of mobile teeth
  - Blocking out of undercuts
  - Adhesive cementation of translucent ceramic and composite restorations
- **Tetric N-Ceram**
  - Class I-V restorations
  - Restoration of deciduous teeth
  - Extended fissure sealing
  - Direct veneers
  - Splinting of mobile teeth
  - Blocking out of undercuts
  - Repair of composite and ceramic veneers

### Contraindications

- **N-Etch**
  - If the stipulated working technique cannot be used, the material should not be applied.
- **Tetric N-Bond**
  - If the patient is allergic to any of the components of Tetric N-Bond or if the stipulated working technique cannot be used.
  - Tetric N-Bond must not be used in combination with purely self-curing composite materials. Always activate dual-curing materials with light.
  - Direct pulp capping with Tetric N-Bond
- **Tetric N-Flow/Tetric N-Ceram**
  - If a dry working field cannot be established or the stipulated application technique cannot be applied.
  - If the patient is known to be allergic to any of the ingredients.

### Side effects

- In rare cases, components of Tetric N-Bond/Tetric N-Flow/Tetric N-Ceram may lead to sensitization. The products should not be used in such cases.
- To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector.

### Interactions

- Glass ionomer-based liners may be dissolved by N-Etch.
- Materials containing eugenol or clove oil may inhibit the polymerization of Tetric N-Bond/Tetric N-Flow/Tetric N-Ceram. Consequently, the application of such materials together with the products must be avoided.
- Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

### Application

#### Direct Restorations

1. Select the shade  
For facilitated access to the treatment field, the lip and cheek retractor OptraGate® can be used.
2. Establish sufficient isolation, preferably with a rubber dam (complete isolation is ensured with, for example, the anatomically shaped OptraDam® Plus).
3. Prepare the cavity according to the requirements of the adhesive technique.
4. Clean the cavity with water spray.
5. Air dry the cavity.
6. Apply a pulp protector (calcium hydroxide) if needed; only cover areas close to the pulp and subsequently apply a pressure-resistant cement (e.g. Vivaglass® Liner).
7. Place a matrix (e.g. the selectively thinned OptraMatrix®) and an interdental wedge, if required.
8. Apply N-Etch on the prepared enamel and then flow the etchant onto the prepared dentin. The etchant should be left to react on the enamel for 15 to 30 seconds, and on the dentin for 10 to 15 seconds. After that, remove all N-Etch with a vigorous water spray for at least 5 seconds. Excess moisture should be removed leaving the dentin surface with a slightly glossy wet appearance (wet-bonding).  
Excess moisture should be removed with either air gun, a dry brush, a foam pellet or other lint-free absorbents.  
**Do not overdry the dentin!**  
The reaction time on unprepared enamel (e.g. fissure sealants) is 30 – 60 seconds.

## Use of the Delivery Forms

Tetric N-Bond is available either in a bottle or the VivaPen.

### VivaPen Delivery Form

- Unscrew the cap of the VivaPen.
- Attach the VivaPen brush cannula. By rotating the front end of the pen, the pre-bent cannula can be brought into the position desired for the application of the adhesive (see picture).

- By activating the click mechanism several times, the VivaPen brush cannula is wetted with Tetric N-Bond. Upon wetting the brush cannula, the yellow colour of Tetric N-Bond becomes visible.

- Apply a thick layer of Tetric N-Bond on the enamel and dentin surfaces, using the enclosed application brush. Brush the material gently into the dentin for at least 10 seconds. Avoid insufficient coverage of the cavity and apply additional material, if required. Remove excess material and the solvent by a gentle stream of air so that the adhesive completely covers the enamel and dentin without pooling.
- Light-cure Tetric N-Bond according to the curing time recommendations below.

A shiny tooth surface prior to the application of the composite shows that all surfaces are completely covered.

- Apply Tetric N-Flow in layers of max. 2 mm or 1.5 mm (Dentin shade) and light-cure according to the curing recommendation indicated below.  
Apply Tetric N-Ceram in layers of max. 2 mm or 1.5 mm (Dentin shade) and contour/adapt the material to the cavity walls using a suitable instrument (e.g. OptraSculpt®).  
The OptraContact® instrument can be used to selectively position proximal contacts.  
Polymerize each layer individually according to the table below. Hold the light emission window as closely as possible to the surface of the restorative material.

Light intensity	Tetric N-Ceram	Tetric N-Flow
≥ 500 mW/cm <sup>2</sup>	20 s	20 s
≥ 1,000 mW/cm <sup>2</sup>	10 s	10 s

Curing Light	Material	Turbo program	High Power program	Soft Power program	Low Power program
LEDition	Tetric N-Bond	–	10 s	–	–
	Tetric N-Flow	–	20 s	–	–
	Tetric N-Ceram	–	20 s	–	–
bluephase C8	Tetric N-Bond	–	–	–	10 s
	Tetric N-Flow	–	15 s	20 s	–
	Tetric N-Ceram	–	15 s	20 s	–
bluephase	Tetric N-Bond	–	–	–	10 s
	Tetric N-Flow	–	10 s	15 s	–
	Tetric N-Ceram	–	10 s	15 s	–
bluephase 20i	Tetric N-Bond	–	–	–	10 s
	Tetric N-Flow	5 s	10 s	15 s	–
	Tetric N-Ceram	5 s	10 s	15 s	–

- Finish the restoration with suitable finishers or fine diamonds. Subsequently, check the occlusion. Polish with silicone polishers (e.g. OptraPol® Next Generation).

### Indirect Restorations

- Prepare the ceramic or composite surface according to the manufacturer's instructions.
- Prepare the tooth in the same way as for direct restorations, steps 1-10.  
**Important note:** Avoid pooling as this may compromise the accuracy of fit of the definitive restoration.
- Subsequently, place the restoration using a light/dual-curing luting material (e.g. Variolink II) and light-cure according to the instructions of the manufacturer.

### Note on the VivaPen

- After each use, the VivaPen needs to be disinfected using disinfectant spray or wipes for hygienic reasons. Additionally, the VivaPen can be protected against contamination with the VivaPen sleeve.
- Discard the used VivaPen brush cannula and replace it with the original cap. A new VivaPen brush cannula must be attached for subsequent applications.
- Do not use oxidizing disinfectants.
- Do not immerse the VivaPen in disinfectant solutions.
- Do not use any other brush cannulas than VivaPen brush cannulas.

### Note on N-Etch Jumbo

#### Refilling of empty N-Etch syringes

Empty 2-gram syringes can be conveniently refilled with the contents of the N-Etch Jumbo syringe up to 15 times.

#### Important

- Although the 2-gram syringe is equipped with a plunger safety stop, it is absolutely necessary to wear gloves, protective goggles and a face mask for reasons of safety.
- Remove the screw-on caps from the 30-gram Jumbo syringe and the empty 2-gram syringe.
- Join the tips of the two syringes with a twisting, clockwise motion. Make sure that the two syringes are firmly joined. Note: Do not strip the thread.
- Refill the small syringe by slowly and steadily pulling the plunger of the small syringe. Simultaneously, exert light pressure on the plunger of the jumbo syringe using the other hand to support the refilling process.
- **Avoid overcharging the small syringe.**
- **Under no circumstances exceed the plunger safety stop of the small syringe – there is a risk of acid burn!**
- Do not point the syringes towards people during refilling.
- Do not fill etchant from the 2-gram syringe back to the Jumbo syringe.
- Disconnect the two syringes by an anticlockwise twist.
- Replace the screw-on cap on the Jumbo syringe.
- Attach an application tip to the 2-gram syringe. The syringe is now ready for use. Apply the etchant as described above.
- Use only the N-Etch application tips delivered with this package. The use of other application tips may have an adverse effect on the viscosity of the etchant.
- In order to ensure traceability, only the enclosed 2-gram syringe should be refilled with the contents of the Jumbo syringe delivered with the same package.

- Please indicate the LOT number on the Jumbo syringe in case of any query.

### Special note

1. If fresh Tetric N-Ceram is applied to polymerized material, it must first be roughened and wetted with Heliobond.
2. Tetric N-Flow/Tetric N-Ceram should have room temperature when applied. Cool temperatures render the material difficult to extrude.
3. Do not use oxidizing disinfectants to disinfect syringes and Cavifils.

### Warning N-Etch

N-Etch contains phosphoric acid and is corrosive. When filling the 2-gram syringe, always wear protective goggles and gloves. Avoid contact with the eyes, mucous membrane and skin (goggles are recommended for both patients and operators). After contact with the eyes, immediately rinse with copious amounts of water for at least 15 minutes and consult a physician/ophthalmologist. After skin contact, immediately rinse with water and soap. After contact with sensitive materials (surfaces, clothing), immediately clean with copious amounts of water.

### Warning Tetric N-Bond

- Tetric N-Bond is an irritant. Avoid contact with skin, mucous membranes and eyes. If the product comes in contact with the skin, rinse with water and soap.
- If Tetric N-Bond comes in contact with the eyes, immediately rinse with copious amounts of water for at least 15 minutes and consult a physician/ophthalmologist.
- In rare cases, contact with the skin may lead to a sensitization to methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

### Warning Tetric N-Flow/Tetric N-Ceram

- Avoid contact of unpolymerized material with the skin/mucous membrane and eyes.
- Unpolymerized material may have a slight irritating effect and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

## Storage and hygiene

- Storage temperature: 2-28 °C (36-82 °F)  
N-Etch recommendation: Refrigerate large stocks at 2-8 °C (36-46 °F).
- Close syringes/Cavifils/bottles/VivaPen immediately after usage. Exposure to light causes premature polymerization.
- Expiry date: see note on Cavifils, syringes, bottles, VivaPen and packaging.
- Do not use any of these products after the indicated date of expiration.
-  If Tetric N-Ceram or Tetric N-Flow is applied from the Cavifil directly in the mouth of the patient, the Cavifil must not be used for more than one patient due to hygienic reasons (prevention of cross-contamination between patients).

## Keep out of the reach of children.

### For use in dentistry only.

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes.

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## Ivoclar Vivadent UK Limited
## Ivoclar Vivadent, Inc.